


ESCMID Global Machine learning workshop 2025: Installing Orange – macOS

Last updated 20/02/2025

This document provides a brief description of how to install the orange data mining software on your Mac computer/laptop. Any problems with installing Orange, please contact me at benjamin.mcfadden@research.uwa.edu.au and I will personally assist you.

1- Visit <https://orangedatamining.com/>



Jan 20, 2025


Foundation models in Orange with just a few lines of code

Foundation models can be used in Orange in just a few lines of code. Python Script can unlock every model from Hugging Face if you are brave enough.


Dec 03, 2024

Partial Least Squares in the Pharmaceutical Industry

Partial Least Squares (PLS) enhances pharmaceutical research by analyzing complex, high-dimensional datasets, addressing multicollinearity, and building predictive models. Orange streamlines its



- 2- After clicking the “Download Orange 3.38.1” button in the previous image, you will see the following screen.



Suggested Download

Benjamin

Orange 3.38.1 for Mac OS

Windows

Standalone installer (default)

⬇️ [Orange3-3.38.1-Miniforge-x86_64.exe](#)

Can be used without administrative privileges.

Portable Orange

⬇️ [Orange3-3.38.1.zip](#)

No installation needed. Just extract the archive and open the shortcut in the extracted folder.

macOS

Orange for Apple silicon

⬇️ [Orange3-3.38.1-Python3.11.8-arm64.dmg](#)

Orange for Intel


⬇️ [Orange3-3.38.1-Python3.10.11-x86_64.dmg](#)

Installing add-ons

Additional features can be added to Orange by installing add-ons. You can find add-on manager in Options menu.

Not sure which installer to select? Click the Apple logo in the top-left corner of your screen, select About This Mac, and check the Chip or Processor field. If you see Apple, select the Orange for Apple Silicon installer. If you see Intel, select the Orange for Intel.

- 3- After clicking the appropriate macOS installation (or the suggested download outlined in red), a disk image will be downloaded.




Name	Size	Kind	Date Added
Orange3-3.38.1-Python3.11.8-arm64.dmg	374.4 MB	Disk Image	Today at 5:58 pm


- 4- When prompted, drag and drop the Orange3 icon to the applications folder.



- 5- Click open when the following appears on screen or when you first open the application.



6- The orange interface should then open.



7- The Orange icon to open the application can also be found on your system.

